


Plastic 12 Channel Fully Autoclavable Micropipette



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic 12 Channel Fully Autoclavable Micropipette Manufactured by Plastic Labware India.

Specification:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

* Lower assembly rotates 360 for optimal working position.

* Optimum piston-spring balance for soft plunger movement.



Plastic 12 Channel Fully Autoclavable Micropipette

Model No:- MPMC 301



- * Easy Volume setting.
- * Easy sequential tip ejection system.
- * Easy user calibration.

We are Manufacturer Exporters and Suppliers of Plastic 12 Channel Fully Autoclavable Micropipette at wholesale competitive prices from India.

<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPMC 301	0.5 - 10 ul	4.00	4.00	1 Pc.
		2.50	2.50	
		1.50	1.50	
MPMC 302	2.0 - 20 ul	4.00	4.00	1 Pc.
		1.50	1.00	
		1.50	1.00	
MPMC 303	5.0 - 50 ul	3.00	2.00	1 Pc.
		1.50	1.00	
		0.50	0.50	
MPMC 304	10 - 100 ul	3.00	1.50	1 Pc.
		1.00	0.50	
		1.00	0.50	
MPMC 305	20 - 200 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	
MPMC 306	30 - 300 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	