

Estd. 1966
D&B
Decide with Confidence

Estd. 1966
DESCO
DELIVERING TO THE WORLD

Plastic 12 Channel Semi Autoclavable Micropipette



www.plasticlabwareindia.com www.descolab.com www.descolaboratoryindia.com

Plastic 12 Channel Semi Autoclavable Micropipette Manufactured by Plastic Labware India.

Specification:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

- * Lower assembly rotates 360 for optimal working position.
- * Optimum piston-spring balance for soft plunger movement.



Plastic 12 Channel Semi Autoclavable Micropipette

Model No:- MPMC 401



- * Easy Volume setting.
- * Easy sequential tip ejection system.
- * Easy user calibration.

We are Manufacturer Exporters and Suppliers of Plastic 12 Channel Semi Autoclavable Micropipette at wholesale competitive prices from India.

<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPMC 401	0.5 - 10 ul	4.0	4.00	1 Pc.
		2.5	2.50	
		1.5	1.50	
MPMC 402	2.0 - 20 ul	4.0	4.00	1 Pc.
		1.5	1.00	
		1.5	1.00	
MPMC 403	5.0 - 50 ul	3.0	2.00	1 Pc.
		1.5	1.00	
		0.5	0.50	
MPMC 404	10 - 100 ul	3.0	1.50	1 Pc.
		1.0	0.50	
		1.0	0.50	
MPMC 405	20 - 200 ul	2.0	0.80	1 Pc.
		1.0	0.50	
		0.7	0.25	
MPMC 406	30 - 300 ul	2.0	0.80	1 Pc.
		1.0	0.50	
		0.7	0.25	