

Estd. 1966
DESCO
DELIVERING TO THE WORLD

D&B
Decide with Confidence

Plastic 8 Channel Fully Autoclavable Micropipette



www.plasticlabwareindia.com www.descolab.com www.descolaboratoryindia.com

Plastic 8 Channel Fully Autoclavable Micropipette Manufactured by Plastic Labware India.

Specification:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

- * Lower assembly rotates 360 for optimal working position.
- * Optimum piston-spring balance for soft plunger movement.



Plastic 8 Channel Fully Autoclavable Micropipette

Model No:- MPMC 101



- * Easy Volume setting.
- * Easy sequential tip ejection system.
- * Easy user calibration.

We are Manufacturer Exporters and Suppliers of Plastic 8 Channel Fully Autoclavable Micropipette at wholesale competitive prices from India.

| <u>Model No.</u> | <u>Volume Range</u> | <u>Accuracy %</u> | <u>Precision (CV%)</u> | <u>Pack Size</u> |
|------------------|---------------------|-------------------|------------------------|------------------|
| MPMC 101 | 0.5 - 10 ul | 4.00 | 4.00 | 1 Pc. |
| | | 2.50 | 2.50 | |
| | | 1.50 | 1.50 | |
| MPMC 102 | 2.0 - 20 ul | 4.00 | 4.00 | 1 Pc. |
| | | 1.50 | 1.00 | |
| | | 1.50 | 1.00 | |
| MPMC 103 | 5.0 - 50 ul | 3.00 | 2.00 | 1 Pc. |
| | | 1.50 | 1.00 | |
| | | 0.50 | 0.50 | |
| MPMC 104 | 10 - 100 ul | 3.00 | 1.50 | 1 Pc. |
| | | 1.00 | 0.50 | |
| | | 1.00 | 0.50 | |
| MPMC 105 | 20 - 200 ul | 2.00 | 0.80 | 1 Pc. |
| | | 1.00 | 0.50 | |
| | | 0.70 | 0.25 | |
| MPMC 106 | 30 - 300 ul | 2.00 | 0.80 | 1 Pc. |
| | | 1.00 | 0.50 | |
| | | 0.70 | 0.25 | |