

Estd. 1966  
**DESCO**  
DELIVERING TO THE WORLD

**D&B**  
Decide with Confidence

## Plastic 8 Channel Semi Autoclavable Micro Pipette



[www.plasticlabwareindia.com](http://www.plasticlabwareindia.com)   [www.descolab.com](http://www.descolab.com)   [www.descolaboratoryindia.com](http://www.descolaboratoryindia.com)

## Plastic 8 Channel Semi Autoclavable Micropipette Manufactured by Plastic Labware India.

### Specification:-

This is Most accurate, reliable, comfortable for extensive pipetting in today's research. Research today demands only highest quality when it comes to pipetting and precision. Handypipette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

- \* Lower assembly rotates 360 for optimal working position.
- \* Optimum piston-spring balance for soft plunger movement.





# Plastic 8 Channel Semi Autoclavable Micropipette

## Model No:- MPMC 201



- \* Easy Volume setting.
- \* Easy sequential tip ejection system.
- \* Easy user calibration.

We are Manufacturer Exporters and Suppliers of Plastic 8 Channel Semi Autoclavable Micropipette at wholesale competitive prices from India.

<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPMC 201	0.5 - 10 ul	4.00	4.00	1 Pc.
		2.50	2.50	
		1.50	1.50	
MPMC 202	2.0 - 20 ul	4.00	4.00	1 Pc.
		1.50	1.00	
		1.50	1.00	
MPMC 203	5.0 - 50 ul	3.00	2.00	1 Pc.
		1.50	1.00	
		0.50	0.50	
MPMC 204	10 - 100 ul	3.00	1.50	1 Pc.
		1.00	0.50	
		1.00	0.50	
MPMC 205	20 - 200 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	
MPMC 206	30 - 300 ul	2.00	0.80	1 Pc.
		1.00	0.50	
		0.70	0.25	