



Decide with Confidence



DELIVERING TO THE WORLD

Plastic Bottle Top Dispenser PP Valve



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic Bottle Top Dispenser PP Valve Manufactured by Plastic Labware India.

Specification:-

Bottle Top Dispenser PP Valve is easy to use, simple construction and smooth, effortless plunger movement for bubble free dispensing.



Plastic Bottle Top Dispenser PP Valve Model No:- BOBD 101



Piston : Specially designed PP piston with silicone O-ring ensures smooth, effortless piston movement and high accuracy.

Each instrument comes with five adaptors to fit in most of the laboratory reagent bottles

standard thread : 30mm

Adaptors : 28mm, 32mm, 36mm, 40mm, and 45mm.

Nozzle cap prevents any unwanted drop on the work spece.

Easy to dis-assemble for cleaning and servicing.

Glass Barrel is protected by a polypropylene sleeve. It also ensures no glass splitting in case of any accidental breakage.

Valve Manifold : Specially designed, single piece Polyproplene valve manifold in the premium model ensures excellent reproducibility.

Instrument is fully autoclavable at 121oC , 15psi. No need for des-assembling the instrument for autoclaving.

Screw type volume adjustment knob with large pointer ensures precise, easy quick volume adjustment.

Calibration : Specially designed universal calibration tool is provided for convenient and quick In-lab user re-calibration. This is in compliance with GLO/ISO Norms.

We are Manufacturer Exporters and Suppliers of Plastic Bottle Top Dispenser PP Valve at wholesale competitive prices from India.

<u>Model No.</u>	<u>Vol Range</u>	<u>Increment</u>	<u>Pack Size</u>
BOBD 101	0.25-2.5 ML	0.05 MI	1 Pc.
BOBD 102	0.5 - 5 ML	0.1 MI	1 Pc.
BOBD 103	1 -10 ML	0.2 MI	1 Pc.
BOBD 104	2.5 -25 ML	0.5 MI	1 Pc.
BOBD 105	5 - 50 ML	1.0 MI	1 Pc.