

Plastic Bottle Top Dispenser PP Valve Model No:- BOBD 101



Decide with Confidence

Plastic Bottle Top Dispenser PP Valve

Estd. 1966



DELIVERING TO THE WORLD



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic Bottle Top Dispenser PP Valve Manufactured by Plastic Labware India.

Specification:-

Bottle Top Dispenser PP Valve is easy to use, simple construction and smooth, effortless plunger movement for bubble free dispensing.



Plastic Bottle Top Dispenser PP Valve Model No:- BOBD 101



Piston : Specially designed PP piston with silicone O-ring ensures smooth, effortless piston movement and high accuracy.

Each instrument comes with five adaptors to fit in most of the laboratory reagent bottles

standard thread : 30mm

Adaptors : 28mm, 32mm, 36mm, 40mm, and 45mm.

Nozzle cap prevents any unwanted drop on the work spece.

Easy to dis-assemble for cleaning and servicing.

Glass Barrel is protected by a polypropylene sleeve. It also ensures no glass splitting in case of any accidental breakage.

Valve Manifold : Specially designed, single piece Polyproplene valve manifold in the premium model ensures excellent reproducibility.

Instrument is fully autoclavable at 121oC , 15psi. No need for des-assembling the instrument for autoclaving.

Screw type volume adjustment knob with large pointer ensures precise, easy quick volume adjustment.

Calibration : Specially designed universal calibration tool is provided for convenient and quick In-lab user re-calibration. This is in compliance with GLO/ISO Norms.

We are Manufacturer Exporters and Suppliers of Plastic Bottle Top Dispenser PP Valve at wholesale competitive prices from India.

<u>Model No.</u>	<u>Vol Range</u>	<u>Increment</u>	<u>Pack Size</u>
BOBD 101	0.25-2.5 ML	0.05 MI	1 Pc.
BOBD 102	0.5 - 5 ML	0.1 MI	1 Pc.
BOBD 103	1 -10 ML	0.2 MI	1 Pc.
BOBD 104	2.5 -25 ML	0.5 MI	1 Pc.
BOBD 105	5 - 50 ML	1.0 MI	1 Pc.