

Estd. 1966  
**DESCO**<sup>®</sup>  
DELIVERING TO THE WORLD

**D&B**  
Decide with Confidence

## Plastic Cell and Tissue Culture Flask



[www.plasticlabwareindia.com](http://www.plasticlabwareindia.com)   [www.descolab.com](http://www.descolab.com)   [www.descolaboratoryindia.com](http://www.descolaboratoryindia.com)

**Plastic Cell and Tissue Culture Flask Manufactured by Plastic Labware India.**

### **Specification:-**

Packaged in sterile, zip-sealable bags, high clarity, 100% virgin polystyrene  
Frosted writing and clear graduations, notched bottom for slip free stacking

Sterilized by E Beam

Clear lot number for batch traceability

Non-pyrogenic, DNase/RNase FREE

We are Manufacturer Exporters and Suppliers of Plastic Cell and Tissue Culture  
Flask at wholesale competitive prices from India.

# Plastic Cell and Tissue Culture Flask Model No:- MPCF 101



<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
MPCF 101	50 ML Sterile Standard	200 Pc.
MPCF 102	50 ML Sterile Vent	200 Pc.
MPCF 103	250 ML Sterile Standard	100 Pc.
MPCF 104	250 ML Sterile Vent	100 Pc.
MPCF 105	750 ML Sterile Standard	40 Pc.
MPCF 106	750 ML Sterile Vent	40 Pc.
MPCF 107	950 ML Sterile Standard	25 Pc.
MPCF 108	950 ML Sterile Vent	25 Pc.