




Plastic Centrifuge Tube Conical Bottom Sterile Model No:- CPCT 201



Decide with Confidence



Plastic Centrifuge Tube Conical Bottom Sterile



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

**Plastic Centrifuge Tube Conical Bottom Sterile Manufactured
by Plastic Labware India.**

Plastic Centrifuge Tube Conical Bottom Sterile Model No:- CPCT 201



Specification:-

Centrifuge Tube are made of clear polypropylene and comes with high density polyethylene caps packed in high quality ziplock bags.

Centrifuge tubes have large writing area for sample identification.

Temperature range from -80°C to 120°C

Leak Proof Design

Easy On/easy Off cap

Clear Accurate graduations

Large white marking spots with reliable quality.

Certified RNase-/DNase-free

We are Manufacturer Exporters and Suppliers of Plastic Centrifuge Tube Conical Bottom Sterile at wholesale competitive prices from India.

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
CPCT 201	15 ML Sterile	500 Pc.
CPCT 202	15 ML Racked	500 Pc.
CPCT 203	15 ML Racked Sterile	500 Pc.
CPCT 204	50 ML Sterile	500 Pc.
CPCT 205	50 ML Racked	300 Pc.
CPCT 206	50 ML Racked Sterile	300 Pc.