



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

Plastic Desiccator Plain



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic Desiccator Plain Manufactured by Plastic Labware India.

Specification:-

Desiccator Plain are made of Polypropylene and Polycarbonate.

These autoclavable plain desiccators can hold vaccum upto 740 mm of Hg for 24 hrs. without any gressing.



Plastic Desiccator Plain Model No:- DEDP 101



The top dome, Moulded in rigid and transparent Polycarbonate, gives a crystal clear view of the desiccant places inside. The Knurled knob on the top provides easy handling of the dome. The lower part is made of Polypropylene. An internal groove is provided on the flange to hold a silicon rubber O-ring. It comes without stop cock. Grease should be applied on the flanges to make it air tight. We are Manufacturer Exporters and Suppliers of Plastic Desiccator Plain at wholesale competitive prices from India.

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
DEDP 101	150 MM	1 Pc.
DEDP 102	200 MM	1 Pc.
DEDP 103	250 MM	1 Pc.
DEDP 104	300 MM	1 Pc.