



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

Plastic Desiccator Vacuum



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic Desiccator Vacuum Manufactured by Plastic Labware India.

Specification:-

Desiccator Vacuum are made of Polypropylene and Polycarbonate.

These autoclavable Vacuum desiccators can hold Vacuum upto 740 mm of Hg for 24 hrs. without any gressing. The top dome, Moulded in rigid and transparent Polycarbonate, gives a crystal clear view of the desiccant places inside.

The Knurled knob on the top provides easy handling of the dome.

The lower part is made of Polypropylene.

Estd. 1966



DELIVERING TO THE WORLD

Plastic Desiccator Vacuum

Model No:- DEDV 101



An internal groove is provided on the flange to hold a silicon rubber O-ring. The Polypropylene stopcock is fitted with a self lubricating PTFE plug. This plug works three way i.e. Vacuum creation, shutting off & Vacuum releasing. We are Manufacturer Exporters and Suppliers of Plastic Desiccator Vacuum at wholesale competitive prices from India.

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
DEDV 101	150 MM	1 Pc.
DEDV 102	200 MM	1 Pc.
DEDV 103	250 MM	1 Pc.
DEDV 104	300 MM	1 Pc.