



# Plastic Egg Shape Magnetic Stirring Bars

## Model No:- MPMB 251



**D&B**  
Decide with Confidence

**DESCO**  
Estd. 1966  
DELIVERING TO THE WORLD

### Plastic EGG Shape Magnetic Stirring Bars

[www.plasticlabwareindia.com](http://www.plasticlabwareindia.com)   [www.descolab.com](http://www.descolab.com)   [www.descolaboratoryindia.com](http://www.descolaboratoryindia.com)

**Plastic Egg Shape Magnetic Stirring Bars Manufactured by Plastic Labware India.**

### Specification:-

These egg-shaped stirring bars are designed especially for use in round-bottom vessels as they spin unrestricted by the vessel contour. Particularly well suited for round bottom flasks; shape ensures complete mixing Offers minimal contact when used in plastic containers Egg shape

We are Manufacturer Exporters and Suppliers of Plastic Egg Shape Magnetic Stirring Bars at wholesale competitive prices from India.

Estd. 1966



®

# Plastic Egg Shape Magnetic Stirring Bars



DELIVERING TO THE WORLD

## Model No:- MPMB 251

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
MPMB 251	5 x 10 MM	50 Pc.
MPMB 252	6 x 15 MM	50 Pc.
MPMB 253	10 x 20 MM	50 Pc.
MPMB 254	12 x 25 MM	25 Pc.
MPMB 255	16 x 30 MM	25 Pc.
MPMB 256	16 x 35 MM	25 Pc.
MPMB 257	20 x 40 MM	10 Pc.
MPMB 258	20 x 50 MM	10 Pc.
MPMB 259	20 x 64 MM	10 Pc.
MPMB 260	20 x 70 MM	10 Pc.