


Plastic Micropipette Mini Fixed Volume Model No:- MPMI 101



Plastic Micropipette Mini Fixed Volume Manufactured by Plastic Labware India.

Specification:-

Micropipette Mini Fixed Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision.

Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

* Light weight and ergonomic design for user comfort.

* Each pipette is individually tested and supplied with a calibration certificate of quality.



Plastic Micropipette Mini Fixed Volume Model No:- MPMI 101



- * Autoclavable tip cone.
- * Mini micropipette is fully autoclavable at 121C & 15 psi.

(Strict Quality Control)

We are Manufacturer Exporters and Suppliers of Plastic Micropipette Mini Fixed Volume at wholesale competitive prices from India.

<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPMI 101	3 ul	1.5	1.0	1 Pc.
MPMI 102	5 ul	1.5	1.0	1 Pc.
MPMI 103	10 ul	1.0	1.0	1 Pc.
MPMI 104	20 ul	0.5	0.5	1 Pc.
MPMI 105	25 ul	0.5	0.5	1 Pc.
MPMI 106	50 ul	0.4	0.3	1 Pc.
MPMI 107	100 ul	0.3	0.3	1 Pc.
MPMI 108	200 ul	0.5	0.5	1 Pc.
MPMI 109	250 ul	0.5	0.5	1 Pc.
MPMI 110	500 ul	0.5	0.5	1 Pc.
MPMI 111	1000 ul	0.5	0.5	1 Pc.