

Plastic Micropipette Non Autoclavable Variable Volume Model No:- MPVA 301



Plastic Micropipette Non Autoclavable Variable Volume Manufactured by Plastic Labware India.

Specification:-

Micropipette Non Autoclavable Variable Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision. Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

- * Individual finely polished stainless steel piston ensure Accuracy
- * Light weight and ergonomic design for user comfort.



Plastic Micropipette Non Autoclavable Variable Volume

Model No:- MPVA 301



- * Quick click volume setting, separate tip ejector and clear digital read out.
- * Each pipette is individually tested and supplied with a calibration certificate of quality.
- * Autoclavable tip cone.
- * Soft touch tip ejection Super blow out for micro volume pipetting (0.5 μ l to 20 μ l models) easy in recalibration and maintenance.

The slender, ribbed hand grip fits comfortably in either right or left hand. We are Manufacturer Exporters and Suppliers of Plastic Micropipette Non Autoclavable Variable Volume at wholesale competitive prices from India.

<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPVA 301	0.2-2.0 ul	12.0	6.0	1 Pc.
		2.5	1.5	
		1.5	0.7	
MPVA 302	0.5 - 10 ul	5.0	2.8	1 Pc.
		1.5	0.8	
		1.0	0.4	
MPVA 303	2.0 - 20 ul	5.0	1.5	1 Pc.
		1.0	0.6	
		1.0	0.3	
MPVA 304	5.0 - 50 ul	2.5	6.0	1 Pc.
		0.7	1.5	
		0.7	0.7	
MPVA 305	10 - 100 ul	2.5	2.8	1 Pc.
		0.8	0.8	
		0.8	0.4	
MPVA 306	20 - 100 ul	2.5	1.5	1 Pc.
		1.0	0.6	
		0.6	0.3	



Plastic Micropipette Non Autoclavable Variable Volume

Model No:- MPVA 301



<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPVA 307	100 - 1000 ul	1.6	0.3	1 Pc.
		0.6	0.2	
		0.5	0.2	
MPVA 308	200 - 2000 ul	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	
MPVA 309	0.5 - 5 ml	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	
MPVA 310	1 - 10 ml	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	