

Plastic Micropipette Semi Autoclavable Variable Volume Model No:- MPVA 201

D&B
Decide with Confidence

Estd. 1966
DESCO
DELIVERING TO THE WORLD

Plastic Micropipette Semi Autoclavable Variable Volume



www.plasticlabwareindia.com www.descolab.com www.descolaboratoryindia.com

Plastic Micropipette Semi Autoclavable Variable Volume Manufactured by Plastic Labware India.

Specification:-

Micropipette Semi Autoclavable Variable Volume is the most accurate, reliable, comfortable for extensive pipetting in today's research.

Research demands only highest quality when it comes to pipetting and precision. Handypette meets these high quality demands by providing exceptionally accurate pipetting while standing up to the rigorous of everyday use.

- * Individual finely polished stainless steel piston ensure Accuracy
- * Light weight and ergonomic design for user comfort.



Plastic Micropipette Semi Autoclavable Variable Volume Model No:- MPVA 201



- * Quick click volume setting, separate tip ejector and clear digital read out.
- * Each pipette is individually tested and supplied with a calibration certificate of quality.
- * Autoclavable tip cone.
- * Soft touch tip ejection Super blow out for micro volume pipetting (0.5 µl to 20 µl models) easy in recalibration and maintenance.

The slender, ribbed hand grip fits comfortably in either right or left hand. We are Manufacturer Exporters and Suppliers of Plastic Micropipette Semi Autoclavable Variable Volume at wholesale competitive prices from India

<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPVA 201	0.2-2.0 ul	12.0	6.0	1 Pc.
		2.5	1.5	
		1.5	0.7	
MPVA 202	0.5 - 10 ul	5.0	2.8	1 Pc.
		1.5	0.8	
		1.0	0.4	
MPVA 203	2.0 - 20 ul	5.0	1.5	1 Pc.
		1.0	0.6	
		1.0	0.3	
MPVA 204	5.0 - 50 ul	2.5	6.0	1 Pc.
		0.7	1.5	
		0.7	0.7	
MPVA 205	10 - 100 ul	2.5	2.8	1 Pc.
		0.8	0.8	
		0.8	0.4	
MPVA 206	20 - 100 ul	2.5	1.5	1 Pc.
		1.0	0.6	
		0.6	0.3	



Plastic Micropipette Semi Autoclavable Variable Volume Model No:- MPVA 201



<u>Model No.</u>	<u>Volume Range</u>	<u>Accuracy %</u>	<u>Precision (CV%)</u>	<u>Pack Size</u>
MPVA 207	100 - 1000 ul	1.6	0.3	1 Pc.
		0.6	0.2	
		0.5	0.2	
MPVA 208	200 - 2000 ul	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	
MPVA 209	0.5 - 5 ml	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	
MPVA 210	1 - 10 ml	1.5	0.3	1 Pc.
		0.6	0.2	
		0.6	0.2	