



Plastic Microtips in Rack Sterile Filter Barrier Tips Manufactured by Plastic Labware India.

Specification:-

Microtips in Rack Sterile is unique boxes are moulded in Polypropylene and are therefore autoclavable.

The Box it self is clear in providing an excellent view of the tips inside while a table like structure is provided inside to hold the tips.

The tip holding spaces are uniformly placed to facilitate the use of multi channel pipettes.

The cover is joined to the base through ball catcher hinges and provides excellent protection to the tips.

Universal Graduation Low Retention Filter Barrier Tip is packed in Racks

We are Manufacturer Exporters and Suppliers of Plastic Microtips in Rack Sterile Filter Barrier Tips at wholesale competitive prices from India



Plastic Microtips in Rack Sterile Filter Barrier Tips Model No:- MPTB 301



<u>Model No.</u>	<u>Capacity</u>	<u>Places</u>	<u>Colour</u>	<u>Pack Size</u>
MPTB 301	10 ul	96 Microtips	White	10 Racks
MPTB 302	20 ul	96 Microtips	White	10 Racks
MPTB 303	50 ul	96 Microtips	White	10 Racks
MPTB 304	100 ul	96 Microtips	White	10 Racks
MPTB 305	200 ul	96 Microtips	White	10 Racks
MPTB 306	1000 ul	100 Microtips	White	6 Racks