

Estd. 1966
DESCO
DELIVERING TO THE WORLD

D&B
Decide with Confidence

Plastic Octagonal Magnetic Stirring Bars

www.plasticlabwareindia.com www.descolab.com www.descolaboratoryindia.com

Plastic Octagonal Magnetic Stirring Bars Manufactured by Plastic Labware India.

Specification:-

Octagonal Magnetic Stirring Bars with Moulded pivot ring angular surface creates high turbulence even at low speed, provides greater area and added turbulence compared to cylindrical bars.

We are Manufacturer Exporters and Suppliers of Plastic Octagonal Magnetic Stirring Bars at wholesale competitive prices from India.

Estd. 1966



®

Plastic Octagonal Magnetic Stirring Bars



DELIVERING TO THE WORLD

Model No:- MPMB 151

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
MPMB 151	8 x 13 MM	50 Pc.
MPMB 152	8 x 16 MM	50 Pc.
MPMB 153	8 X 19 MM	50 Pc.
MPMB 154	8 x 22 MM	50 Pc.
MPMB 155	8 x 25 MM	50 Pc.
MPMB 156	8 x 30 MM	50 Pc.
MPMB 157	8 x 38 MM	50 Pc.
MPMB 158	8 x 41 MM	50 Pc.
MPMB 159	8 x 51 MM	50 Pc.
MPMB 160	8 x 64 MM	50 Pc.
MPMB 161	9.5 x 13 MM	25 Pc.
MPMB 162	9.5 x 19 MM	25 Pc.
MPMB 163	9.5 x 25 MM	25 Pc.
MPMB 164	9.5 x 38 MM	25 Pc.
MPMB 165	9.5 x 41 MM	25 Pc.
MPMB 166	9.5 x 51 MM	25 Pc.
MPMB 167	13 x 50 MM	10 Pc.
MPMB 168	13 x 75 MM	10 Pc.