



Plastic Petri Dish Culture One Compartment Model No:- PD1C 101



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

Plastic Petri Dish Culture One Compartment



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

**Plastic Petri Dish Culture One Compartment Manufactured by
Plastic Labware India.**

Specification:-

Petri Dish Culture are transparent and ready to use.

These Culture Petri Dishes are made of optically clear, Non-toxic, Polystyrene, Polypropylene and Polycarbonate

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)



Plastic Petri Dish Culture One Compartment Model No:- PD1C 101



These Petri Dishes are packed individually in polyethylene bags sterilized by radiation.

We are Manufacturer Exporters and Suppliers of Plastic Petri Dish Culture One Compartment at wholesale competitive prices from India.

<u>Model No.</u>	<u>Size & Material</u>	<u>Pack Size</u>
PD1C 101	55 x 15 MM Polystyrene	20 Pairs
PD1C 102	55 x 15 MM Polypropylene	20 Pairs
PD1C 111	85 x 15 MM Polystyrene	20 Pairs
PD1C 112	85 x 15 MM Polypropylene	20 Pairs
PD1C 113	85 x 15 MM Polycarbonate	10 Pairs
PD1C 121	100 x 15 MM Polystyrene	20 Pairs
PD1C 122	100 x 15 MM Polypropylene	20 Pairs
PD1C 123	100 x 15 MM Polycarbonate	10 Pairs
PD1C 131	150 x 15 MM Polystyrene	20 Pairs
PD1C 132	150 x 15 MM Polypropylene	20 Pairs
PD1C 133	150 x 15 MM Polycarbonate	10 Pairs