



Plastic Pipette Controller

Model No:- PPPC 101



Decide with Confidence

Plastic Pipette Controller

Estd. 1966



DELIVERING TO THE WORLD



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic Pipette Controller Manufactured by Plastic Labware India.

Specification:-

- Can fill 25 ml pipette under 4 seconds
- Safety valve and hydrophobic filters provide double protection against fluid penetration

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)



Plastic Pipette Controller

Model No:- PPPC 101



- Two different speed modes-High & Low along with Gravity drain
- The pump speed can also be fine-tuned by varying finger pressure on operating knobs for better control of speed
- Cadmium free environment friendly NiMH batteries
- Batteries can be changed very easily by the user
- The intelligent charger prevents over charging / heating of batteries
- Low battery indicator
- The specially designed Desktop stand as wall mounted stand enables the Pipette Controller to be charged while resting on it.
- The Pipette Controller can hold the serological pipette while resting on Desktop Stand

Technical Data

Power 2 x 1.2 V/ 1000 mAh NiMH batteries

Battery Service Life Eight hours continuous use

Battery Charging Time 2 - 3 hours to charge fully

Volumetric Range Plastic or Glass pipette from 0.1 ml up to 200 ml

Safety System Replaceable hydrophobic 0.45 µm PTFE filter

Safety check valve

Standard Accessory Battery charger

Desktop Stand

Hydrophobic PTFE membrane filter 0.45 µm 2Pcs

Hydrophobic PTFE membrane filter 0.45 µm 1Pcs

Wall mounting stand

We are Manufacturer Exporters and Suppliers of Plastic Pipette Controller at wholesale competitive prices from India.

Model No.

Size

Pack Size

PPPC 101

25 ML Pipette

1 Pc.