

Plastic Scoop

Model No:- ACSC 101



 

Plastic Scoop



www.plasticlabwareindia.com www.descolab.com www.descolaboratoryindia.com

Plastic Scoop Manufactured by Plastic Labware India.

Specification:-

Scoop are excellent replacement for metallic ones.

These are Moulded in polypropylene which are autoclavable and completely eliminate the danger of rusting or corrosion.

We are Manufacturer Exporters and Suppliers of Plastic Scoop at wholesale competitive prices from India.

Estd. 1966



Plastic Scoop

Model No:- ACSC 101

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
ACSC 101	2 ML	12 Pc.
ACSC 102	5 ML	12 Pc.
ACSC 103	10 ML	12 Pc.
ACSC 104	25 ML	12 Pc.
ACSC 105	50 ML	12 Pc.
ACSC 106	100 ML	12 Pc.
ACSC 107	250 ML	12 Pc.
ACSC 108	500 ML	6 Pc.
ACSC 109	1000 ML	6 Pc.
ACSC 110	Scoop- Set/9	1 Set