

Plastic Test Tube Stand 3 Tier PC Model No:- TTTS 501




Decide with Confidence

Plastic Test Tube Stand 3 Tier PC

Estd. 1966



DELIVERING TO THE WORLD



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

**Plastic Test Tube Stand 3 Tier Polycarbonate Manufactured by
Plastic Labware India.**

Plastic Test Tube Stand 3 Tier PC Model No:- TTTS 501



Specification:-

Test Tube Stand 3 Tier PC is made of Polycarbonate ingenious racks with a 3-tier design which provides clear view of the tube contents & convenient insertion & removal of tubes. The PC rack has opaque upper & middle plates.

These autoclavable Racks can with stand sub-freezing temperatures.

We are Manufacturer Exporters and Suppliers of Plastic Test Tube Stand 3 Tier Polycarbonate at wholesale competitive prices from India.

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
TTTS 501	13 MM x 31 Tubes	2 Pc.
TTTS 502	16 MM x 31 Tubes	2 Pc.
TTTS 503	20 MM x 20 Tubes	2 Pc.
TTTS 504	25 MM x 12 Tubes	2 Pc.
TTTS 505	13 MM x 62 Tubes	2 Pc.
TTTS 506	16 MM x 62 Tubes	2 Pc.
TTTS 507	20 MM x 40 Tubes	2 Pc.
TTTS 508	25 MM x 24 Tubes	2 Pc.
TTTS 509	25 MM x 18 Tubes	2 Pc.
TTTS 510	25 MM x 36 Tubes	2 Pc.
TTTS 511	32 MM x 12 Tubes	2 Pc.