

# Plastic Test Tube Stand 3 Tier PP Model No:- TTTS 601




Decide with Confidence

## Plastic Test Tube Stand 3 Tier PP

Estd. 1966



DELIVERING TO THE WORLD



[www.plasticlabwareindia.com](http://www.plasticlabwareindia.com)

[www.descolab.com](http://www.descolab.com)

[www.descolaboratoryindia.com](http://www.descolaboratoryindia.com)

**Plastic Test Tube Stand 3 Tier Polypropylene Manufactured by  
Plastic Labware India.**

# Plastic Test Tube Stand 3 Tier PP Model No:- TTTS 601



## Specification:-

Test Tube Stand 3 Tier PP is made of Polypropylene ingenious racks with a 3-tier design provide clear view of the tube contents & convenient insertion & removal of tubes. The PP rack has opaque upper & middle plates.

These autoclavable Racks can withstand sub-freezing temperatures.

We are Manufacturer Exporters and Suppliers of Plastic Test Tube Stand 3 Tier Polypropylene at wholesale competitive prices from India..

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
TTTS 601	13 MM x 31 Tubes	2 Pc.
TTTS 602	16 MM x 31 Tubes	2 Pc.
TTTS 603	20 MM x 20 Tubes	2 Pc.
TTTS 604	25 MM x 12 Tubes	2 Pc.
TTTS 605	25 MM x 18 Tubes	2 Pc.
TTTS 606	32 MM x 12 Tubes	2 Pc.
TTTS 607	25 MM x 36 Tubes	2 Pc.
TTTS 608	13 MM x 18 Tubes	4 Pc.
TTTS 609	15 MM x 18 Tubes	4 Pc.
TTTS 610	18 MM x 12 Tubes	4 Pc.
TTTS 611	13 MM x 62 Tubes	2 Pc.
TTTS 612	16 MM x 62 Tubes	2 Pc.
TTTS 613	20 MM x 40 Tubes	2 Pc.
TTTS 614	13 MM x 48 Tubes	2 Pc.
TTTS 615	13 MM x 72 Tubes	2 Pc.