



Plastic Test Tube Stand Wire Pattern Model No:- TTTS 301




Decide with Confidence

Plastic Test Tube Stand Wire Pattern

Estd. 1966



DELIVERING TO THE WORLD



www.plasticlabwareindia.com

www.descolab.com

www.descolaboratoryindia.com

Plastic Test Tube Stand Wire Pattern Manufactured by Plastic Labware India.



Plastic Test Tube Stand Wire Pattern Model No:- TTTS 301



Specification:-

Test Tube Stand Wire Pattern is an economical substitute for wire racks.

These submersible autoclavable racks are stackable when empty.

These are available in different colours and are made of a special blend of Polypropylene that makes it sink in a water bath and maintain their stability even in agitated water. These space saving racks can easily be disassembled for cleaning purposes. The base is in the form of a grid with square openings and has fold up sides and end bars that lock into rack top.

We are Manufacturer Exporters and Suppliers of Plastic Test Tube Stand Wire Pattern at wholesale competitive prices from India.

<u>Model No.</u>	<u>Size</u>	<u>Pack Size</u>
TTTS 301	13 MM x 90 Tubes	6 Pc.
TTTS 302	16 MM x 60 Tubes	6 Pc.
TTTS 303	20 MM x 40 Tubes	6 Pc.
TTTS 304	25 MM x 24 Tubes	6 Pc.
TTTS 305	30 MM x 21 Tubes	6 Pc.